

Product Datasheet

Recombinant Human Tubulin-specific chaperone cofactor E-like protein(TBCEL) (orb1652247)

Description	Recombinant Human Tubulin-specific chaperone cofactor E-like protein(TBCEL)
Species/Host	E. coli
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-SUMO-tagged
Note	For research use only
Expiration Date	6 months from date of receipt.