

Product Datasheet

MS4A14 antibody (orb165243)

Description	Rabbit polyclonal antibody to MS4A14 (C-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	MS4A14
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	76580
Uniprot ID	Q96JA4
NCBI	NP_115986.3 , NP_001073160.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

