

Product Datasheet

PRRT4 antibody (orb165248)

Description	Rabbit polyclonal antibody to PRRT4 (C-term).
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	PRRT4
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	92712
Uniprot ID	C9JH25
NCBI	NP_001108198.2 , NP_001167635.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

