

Product Datasheet

Recombinant Escherichia coli Plasmid-derived single-stranded DNA-binding protein(ssbF) (orb1652511)

Description	Recombinant Escherichia coli Plasmid-derived single-stranded DNA-binding protein(ssbF)
Species/Host	E. coli
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-SUMO-tagged
Note	For research use only
Expiration Date	6 months from date of receipt.