

Product Datasheet

DYTN antibody (orb165259)

Description	Rabbit polyclonal antibody to DYTN (C-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	DYTN
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	65320
Uniprot ID	A2CJ06
NCBI	NP_001087199.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

