

Product Datasheet

PLEKHA1 antibody (orb165284)

Description	Rabbit polyclonal antibody to PLEKHA1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	PLEKHA1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	45553
Uniprot ID	Q9HB21
NCBI	NP_001001974.1 , NP_067635.2 , NP_001182537.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

