

Product Datasheet

Recombinant Human Cyclin-dependent kinase 2-associated protein 1(CDK2AP1) (orb1653203)

Description	Recombinant Human Cyclin-dependent kinase 2-associated protein 1(CDK2AP1)
Species/Host	E. coli
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal GST-tagged
Note	For research use only
Expiration Date	6 months from date of receipt.