

Product Datasheet

HAUS8 antibody (orb165323)

Description	Rabbit polyclonal antibody to HAUS8 (N-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	HAUS8
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	44857
Uniprot ID	Q9BT25
NCBI	NP_219485.1 , NP_001011699.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

