

Product Datasheet

Goat Insulin-like growth factor 1 (IGF-1)ELISA Kit (orb1653389)

Description	Goat Insulin-like growth factor 1 (IGF-1)ELISA Kit
Reactivity	Goat
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.625 ng/mL-40 ng/mL. Sensitivity: 0.3 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.