

## Product Datasheet

# Human ACE2 / Angiotensin-converting enzyme 2 ELISA Kit (orb1653394)

|                            |                                                                                |
|----------------------------|--------------------------------------------------------------------------------|
| <b>Description</b>         | Human ACE2 / Angiotensin-converting enzyme 2 ELISA Kit                         |
| <b>Reactivity</b>          | Human                                                                          |
| <b>Tested Applications</b> | ELISA                                                                          |
| <b>Note</b>                | For research use only                                                          |
| <b>Application notes</b>   | Detect Range: 0.78-50 ng/mL.<br>Sensitivity: 0.34 ng/mL. Assay Type: Sandwich. |
| <b>Expiration Date</b>     | Please enquire.                                                                |