

## Product Datasheet

### Human EPRS / Bifunctional glutamate/proline--tRNA ligase ELISA Kit (orb1653402)

|                            |                                                                                  |
|----------------------------|----------------------------------------------------------------------------------|
| <b>Description</b>         | Human EPRS / Bifunctional glutamate/proline--tRNA ligase ELISA Kit               |
| <b>Reactivity</b>          | Human                                                                            |
| <b>Tested Applications</b> | ELISA                                                                            |
| <b>Note</b>                | For research use only                                                            |
| <b>Application notes</b>   | Detect Range: 0.156-10 ng/mL.<br>Sensitivity: 0.055 ng/mL. Assay Type: Sandwich. |
| <b>Expiration Date</b>     | Please enquire.                                                                  |