

## Product Datasheet

# Mouse Lingo1 / Leucine-rich repeat and immunoglobulin-like domain-containing nogo receptor-interacting protein 1 ELISA Kit (orb1653419)

|                            |                                                                                                         |
|----------------------------|---------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Mouse Lingo1 / Leucine-rich repeat and immunoglobulin-like domain-containing nogo receptor-interacti... |
| <b>Reactivity</b>          | Mouse                                                                                                   |
| <b>Tested Applications</b> | ELISA                                                                                                   |
| <b>Note</b>                | For research use only                                                                                   |
| <b>Application notes</b>   | Detect Range: 0.312-20 ng/mL.<br>Sensitivity: 0.11 ng/mL. Assay Type: Sandwich.                         |
| <b>Expiration Date</b>     | Please enquire.                                                                                         |