

Product Datasheet

Human Cation-dependent mannose-6-phosphate receptor (M6PR) ELISA Kit (orb1653511)

Description	Human Cation-dependent mannose-6-phosphate receptor (M6PR) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 1.56-100 ng/mL. Sensitivity: 0.63 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.