

## Product Datasheet

# Human Methionine adenosyltransferase 2 subunit beta (MAT2B) ELISA Kit (orb1653631)

|                            |                                                                                                                                                                                                       |
|----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Human Methionine adenosyltransferase 2 subunit beta (MAT2B) ELISA Kit                                                                                                                                 |
| <b>Reactivity</b>          | Human                                                                                                                                                                                                 |
| <b>Tested Applications</b> | ELISA                                                                                                                                                                                                 |
| <b>Note</b>                | For research use only                                                                                                                                                                                 |
| <b>Application notes</b>   | Detect Range: Request Information. Sensitivity: Request Information. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids. |
| <b>Expiration Date</b>     | Please enquire.                                                                                                                                                                                       |