

Product Datasheet

CAMSAP3 antibody (orb165365)

Description	Rabbit polyclonal antibody to CAMSAP3 (N-term).
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	CAMSAP3
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	134750
Uniprot ID	Q9P1Y5
NCBI	NP_065953.1 , NP_001073898.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

