

Product Datasheet

Cat E3 ubiquitin-protein ligase Mdm2 (MDM2) ELISA Kit (orb1653740)

Description	Cat E3 ubiquitin-protein ligase Mdm2 (MDM2) ELISA Kit
Reactivity	Feline
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: Request Information. Sensitivity: Request Information. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.