

Product Datasheet

Human Putative homeobox protein Meis3-like 2 (MEIS3P2) ELISA Kit (orb1653789)

Description Human Putative homeobox protein Meis3-like 2 (MEIS3P2) ELISA Kit

Reactivity Human

Tested Applications ELISA

Note For research use only

Application notes Detect Range: Request Information. Sensitivity: Request Information. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids.

Expiration Date Please enquire.