

Product Datasheet

GLT8D1 antibody (orb165385)

Description	Rabbit polyclonal antibody to GLT8D1 (C-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	GLT8D1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	41935
Uniprot ID	Q68CQ7
NCBI	NP_001265210.1 , NP_001265209.1 , NP_060916.1 , NP_001010983.1 , NP_690909.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

