

Product Datasheet

UCK1 antibody (orb165394)

Description	Rabbit polyclonal antibody to UCK1 (C-term).
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	UCK1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	31435
Uniprot ID	Q9HA47
NCBI	NP_001248379.1 , NP_001248380.1 , NP_001129426.1 , NP_113620.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

