

Product Datasheet

Cat NADH-ubiquinone oxidoreductase chain 5 (MT-ND5) ELISA Kit (orb1654336)

Description Cat NADH-ubiquinone
oxidoreductase chain 5 (MT-ND5)
ELISA Kit

Reactivity Feline

**Tested
Applications** ELISA

Note For research use only

Application notes Detect Range: Request
Information. Sensitivity: Request
Information. Assay Time: 1-4.5 h.
Detection wavelength: 450 nm.
Sample Volume: 1-200 µL.
Sample Type: Serum, Plasma,
Other biological fluids.

Expiration Date Please enquire.