

Product Datasheet

ZNF33A antibody (orb165453)

Description	Rabbit polyclonal antibody to ZNF33A (N-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	ZNF33A
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	94384
Uniprot ID	Q06730
NCBI	NP_001265104.1 , NP_001265100.1 , NP_001265105.1 , NP_001265099.1 , NP_001265106.1 , NP_001265102.1 , NP_008885.1 , NP_008905.1 , NP_001265103.1 , NP_001265107.1 , NP_001265108.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

