

Product Datasheet

OR2T5 antibody (orb165456)

Description	Rabbit polyclonal antibody to OR2T5 (N-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	OR2T5
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	35595
Uniprot ID	Q6IEZ7
NCBI	NP_001004697.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

