

Product Datasheet

POMZP3 antibody (orb165457)

Description	Rabbit polyclonal antibody to POMZP3 (N-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	POMZP3
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	20620
Uniprot ID	Q6PJE2
NCBI	NP_694537.1 , NP_036362.3
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

