

Product Datasheet

FAHD2B antibody (orb165473)

Description	Rabbit polyclonal antibody to FAHD2B (C-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	FAHD2B
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	34613
Uniprot ID	Q6P2I3
NCBI	NP_955368.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

