

Product Datasheet

Human Polyadenylate-binding protein 2 (PABPN1) ELISA Kit (orb1654814)

Description	Human Polyadenylate-binding protein 2 (PABPN1) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 15.6-1000 pg/mL. Sensitivity: 7.2 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.