

Product Datasheet

Human Protein-arginine deiminase type-3 (PADI3) ELISA Kit (orb1654816)

Description Human Protein-arginine
deiminase type-3 (PADI3) ELISA
Kit

Reactivity Human

**Tested
Applications** ELISA

Note For research use only

Application notes Detect Range: 0.156-10 ng/mL.
Sensitivity: 0.078 ng/mL. Assay
Time: 1-4.5 h. Detection
wavelength: 450 nm. Sample
Volume: 1-200 µL. Sample Type:
Serum, Plasma, Other biological
fluids.

Expiration Date Please enquire.