

## Product Datasheet

### Human Platelet activating factor (PAF) ELISA Kit (orb1654820)

|                            |                                                                                                                                                                                                        |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Human Platelet activating factor<br>(PAF) ELISA Kit                                                                                                                                                    |
| <b>Reactivity</b>          | Human                                                                                                                                                                                                  |
| <b>Tested Applications</b> | ELISA                                                                                                                                                                                                  |
| <b>Note</b>                | For research use only                                                                                                                                                                                  |
| <b>Application notes</b>   | Detect Range: 4.69-300 ng/mL.<br>Sensitivity: 1.17 ng/mL. Assay Time:<br>1-4.5 h. Detection wavelength: 450<br>nm. Sample Volume: 1-200 µL.<br>Sample Type: Serum, Plasma, Other<br>biological fluids. |
| <b>Expiration Date</b>     | Please enquire.                                                                                                                                                                                        |