

Product Datasheet

Bovine Pregnancy Specific Protein B (PSPB) ELISA Kit (orb1654823)

Description	Bovine Pregnancy Specific Protein B (PSPB) ELISA Kit
Reactivity	Bovine
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.06-24 ng/mL. Sensitivity: 0.06 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.