

## Product Datasheet

# Rat Platelet-Derived Growth Factor AB (PDGF-AB) ELISA Kit (orb1654883)

|                            |                                                                                                                                                                                        |
|----------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rat Platelet-Derived Growth Factor AB (PDGF-AB) ELISA Kit                                                                                                                              |
| <b>Reactivity</b>          | Rat                                                                                                                                                                                    |
| <b>Tested Applications</b> | ELISA                                                                                                                                                                                  |
| <b>Note</b>                | For research use only                                                                                                                                                                  |
| <b>Application notes</b>   | Detect Range: 7.8-500 pg/mg. Sensitivity: 1.95 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids. |
| <b>Expiration Date</b>     | Please enquire.                                                                                                                                                                        |