

Product Datasheet

Human Prostaglandin E2 (PG-E2) ELISA Kit (orb1654937)

Description	Human Prostaglandin E2 (PG-E2) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 1.5-300 pg/mL. Sensitivity: 0.75 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.