

Product Datasheet

IL34 antibody (orb165497)

Description	Rabbit polyclonal antibody to IL34 (N-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	IL34
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	27482
Uniprot ID	Q6ZMJ4
NCBI	NP_689669.2 , NP_001166243.1 , NP_001166242.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

