

Product Datasheet

Monkey Progesterone (PROG) ELISA Kit (orb1655333)

Description	Monkey Progesterone (PROG) ELISA Kit
Reactivity	Monkey
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.4-30 ng/mL. Sensitivity: 0.2 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.