

Product Datasheet

FRMD3 antibody (orb165541)

Description	Rabbit polyclonal antibody to FRMD3
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	FRMD3
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	68772
Uniprot ID	A2A2Y4
NCBI	NP_001231888.1 , NP_001231890.1 , NP_777598.3 , NP_001231891.1 , NP_001231889.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

