

Product Datasheet

Human Receptor-type tyrosine-protein phosphatase zeta (PTPRZ1) ELISA Kit (orb1655534)

Description	Human Receptor-type tyrosine-protein phosphatase zeta (PTPRZ1) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 93.75-6000 pg/mL. Sensitivity: 23.44 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.