

Product Datasheet

PDE4C antibody (orb165563)

Description	Rabbit polyclonal antibody to PDE4C (N-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	PDE4C
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	79902
Uniprot ID	Q08493
NCBI	NP_000914.2 , NP_001092288.1 , NP_001092289.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

