

Product Datasheet

Horse Angiotensinogen (AGT) ELISA Kit (orb1655693)

Description Horse Angiotensinogen (AGT) ELISA Kit

Reactivity Equine

Tested Applications ELISA

Note For research use only

Application notes Detect Range: 31.25-2000 pg/mL.
Sensitivity: 15.6 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm.
Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids.

Expiration Date Please enquire.