

Product Datasheet

Rat Agouti Related Protein (AGRP) ELISA Kit (orb1655694)

Description	Rat Agouti Related Protein (AGRP) ELISA Kit
Reactivity	Rat
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 9.38-600 pg/mL. Sensitivity: 3.5 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.