

Product Datasheet

Human Agmatine Ureohydrolase/Agmatinase (AGMAT) ELISA Kit (orb1655700)

Description	Human Agmatine Ureohydrolase/Agmatinase (AGMAT) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 39-2500 pg/mL. Sensitivity: 9.75 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.