

Product Datasheet

Mouse Angiogenic factor with G patch and FHA domains 1 (AGGF1) ELISA Kit (orb1655702)

Description Mouse Angiogenic factor with G patch and FHA domains 1 (AGGF1) ELISA Kit

Reactivity Mouse

Tested Applications ELISA

Note For research use only

Application notes Detect Range: 78.1-5000 pg/mL. Sensitivity: 34 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids.

Expiration Date Please enquire.