

Product Datasheet

Rat Asymmetrical dimethylarginine (ADMA) ELISA Kit (orb1655743)

Description	Rat Asymmetrical dimethylarginine (ADMA) ELISA Kit
Reactivity	Rat
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 12.35-1000 ng/mL. Sensitivity: 4.51 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.