

Product Datasheet

IL17RC antibody (orb165577)

Description	Rabbit polyclonal antibody to IL17RC (N-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	IL17RC
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	86240
Uniprot ID	Q8NAC3
NCBI	NP_001190192.1 , NP_703191.1 , NP_116121.2 , NP_001190194.1 , NP_001190193.1 , NP_703190.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

