

Product Datasheet

Mouse Activin B (ACV-B) ELISA Kit (orb1655795)

Description	Mouse Activin B (ACV-B) ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 6.25-400 pg/mL. Sensitivity: 1.56 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.