



Product Datasheet

RAB5B antibody (orb165580)

Description	Rabbit polyclonal antibody to RAB5B
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	RAB5B
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	23707
Uniprot ID	P61020
NCBI	NP_002859.1 , NP_001238965.1 , NP_001238966.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse liver tis...