

Product Datasheet

KIZUNA antibody (orb165583)

Description	Rabbit polyclonal antibody to KIZUNA (N-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	KIZ
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	75111
Uniprot ID	Q2M2Z5
NCBI	NP_001156494.1 , NP_001263318.1 , NP_060944.3 , NP_001156495.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

