

Product Datasheet

Rat Angiotensin converting enzyme 2 (ACE2) ELISA Kit (orb1655836)

Description	Rat Angiotensin converting enzyme 2 (ACE2) ELISA Kit
Reactivity	Rat
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 1.56-100 IU/mL. Sensitivity: 0.39 IU/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.