

Product Datasheet

Human Tyrosine-protein kinase ABL1 (ABL1) ELISA Kit (orb1655853)

Description Human Tyrosine-protein kinase
ABL1 (ABL1) ELISA Kit

Reactivity Human

Tested ELISA

Applications

Note For research use only

Application notes Detect Range: 25-1600 pg/mL.
Sensitivity: 6.25 pg/mL. Assay
Time: 1-4.5 h. Detection
wavelength: 450 nm. Sample
Volume: 1-200 µL. Sample Type:
Serum, Plasma, Other biological
fluids.

Expiration Date Please enquire.