

Product Datasheet

HIST1H2AG antibody (orb165588)

Description	Rabbit polyclonal antibody to HIST1H2AG (C-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	H2AC11
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	14091
Uniprot ID	P0C0S8
NCBI	NP_066408.1 , NP_003501.1 , NP_003500.1 , NP_003502.1 , NP_003505.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

