

Product Datasheet

Human 12-hydroxyeicosatetraenoic acid (12-HETE) ELISA Kit (orb1655916)

Description	Human 12-hydroxyeicosatetraenoic acid (12-HETE) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 12.35-1000 ng/mL. Sensitivity: 4.64 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.