

Product Datasheet

Mouse Retinol-binding protein 4 (RBP4) ELISA Kit (orb1655926)

Description Mouse Retinol-binding protein 4
(RBP4) ELISA Kit

Reactivity Mouse

Tested ELISA

Applications

Note For research use only

Application notes Detect Range: 0.312-20 ng/mL.
Sensitivity: 0.142 ng/mL. Assay
Time: 1-4.5 h. Detection
wavelength: 450 nm. Sample
Volume: 1-200 µL. Sample Type:
Serum, Plasma, Other biological
fluids.

Expiration Date Please enquire.