

Product Datasheet

PGLYRP1 antibody (orb165609)

Description	Rabbit polyclonal antibody to PGLYRP1 (C-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	PGLYRP1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	21731
Uniprot ID	O75594
NCBI	NP_005082.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

